

## Image Manipulation using Nibabel:

This Jupyter notebook provides a detailed introduction to working with neuroimaging data using the **Nibabel** library in Python. Below are the key points covered in the notebook:

### Overview

The notebook focuses on:

- Loading, modifying, saving, and visualizing neuroimages.
- Understanding the data structures involved in neuroimaging.
- Using **Nibabel**, a low-level Python library for handling various neuroimaging formats (e.g., NIfTI, FreeSurfer, DICOM).

### Key Sections

#### 1. Setup

- **Imports:** Essential libraries such as nilearn, nibabel, numpy, and matplotlib are imported to handle image processing and visualization.

```
python
```

```
from nilearn import plotting
```

```
import numpy as np
```

```
import nibabel as nb
```

```
import pylab as plt
```

```
%matplotlib inline
```

#### 2. Loading and Inspecting Images

- **Loading an Image:** A functional MRI image is loaded using nibabel.load(). The image file is in NIfTI format (.nii.gz).

```
python
```

```
img = nb.load('/data/ds000114/sub-01/ses-test/func/sub-01_ses-test_task-fingerfootlips_bold.nii.gz')
```

- **Inspecting Image Metadata:** The header and affine transformation matrix of the image are printed. These provide details about the image's dimensions, voxel sizes, and orientation.

python

```
print(img) # Prints metadata such as data shape, affine matrix, and header information.
```

### 3. Accessing Specific Parameters

- **Extracting Data:** The image data is extracted into a NumPy array using `img.get_fdata()`. This array can be manipulated like any standard NumPy array.

python

```
data = img.get_fdata()
```

```
affine = img.affine # Affine transformation matrix
```

```
header = img.header['pixdim'] # Voxel resolution and TR (repetition time)
```

### 4. Visualizing Data

- **Plotting Slices:** A slice of the data is plotted using `matplotlib.pyplot.imshow()`. The shape of the data is also printed.

python

```
plt.imshow(data[:, :, data.shape[2] // 2, 0].T, cmap='Greys_r')
```

```
print(data.shape) # (64, 64, 30, 184)
```

### 5. Exercise: Loading T1 Data

- An exercise is provided to load a T1-weighted anatomical image and plot it similarly to the functional image.

python

```
t1 = nb.load('/data/ds000114/sub-01/ses-test/anat/sub-01_ses-test_T1w.nii.gz')
```

```
data = t1.get_fdata()
```

```
plt.imshow(data[:, :, data.shape[2] // 2].T, cmap='Greys_r')
```

```
print(data.shape) # (256, 156, 256)
```

### 6. Affine Transformation

- The affine matrix is used to convert voxel coordinates into real-world coordinates (in millimeters). Additionally, it encodes information about voxel sizes and axis orientation.

python

```
x, y, z, _ = np.linalg.pinv(affine).dot(np.array([0, 0, 0, 1])).astype(int)

nb.aff2axcodes(affine) # Returns axis orientation ('L', 'A', 'S')

nb.affines.voxel_sizes(affine) # Returns voxel sizes ([3.99999995, 4.00000009, 3.99997491])
```

## 7. Header Information

- The header contains metadata about the image (e.g., voxel sizes, units). It can be accessed directly or through helper functions.

python

```
t1.header.get_zooms() # Returns voxel size in each dimension (1.0, 1.2993759, 1.0)

t1.header.get_xyz_t_units() # Returns units ('mm', 'sec')

t1.header.get_qform() # Returns qform matrix for spatial orientation
```

## 8. Creating and Saving Images

- The notebook demonstrates how to modify image data (e.g., rescaling) and save it as a new NIfTI file while preserving the original header and affine.

python

```
# Rescale data to unsigned byte format (uint8)

rescaled = ((data - data.min()) * 255. / (data.max() - data.min())).astype(np.uint8)

# Save new image with same affine and header

new_img = nb.Nifti1Image(rescaled, affine=img.affine, header=img.header)

nb.save(new_img, '/tmp/rescaled_image.nii.gz')
```

- The notebook also shows how to correct issues with the header when saving images by using `set_data_dtype()`.

python

```
img.set_data_dtype(np.uint8) # Corrects header dtype before saving again.
```

## 9. File Size Comparison

- After rescaling the image data type from int16 to uint8, the file size is reduced significantly.

bash

```
!du -hL /tmp/rescaled_image.nii.gz /data/ds000114/sub-01/ses-test/func/sub-01_ses-test_task-fingerfootlips_bold.nii.gz
```

*# Output shows that the rescaled file is smaller.*

## **Conclusion**

The notebook introduces key functionalities of Nibabel for working with neuroimaging data:

- Loading and inspecting NIfTI images.
- Accessing metadata such as headers and affine transformations.
- Visualizing slices of neuroimaging data.
- Modifying and saving images while preserving important metadata like headers and affines.

This foundational knowledge prepares users for more advanced neuroimaging analysis workflows using libraries like Nibabel and Nilearn.